





Single-use



Prone position headrest with mirror

For pressure-relieving positioning of patients



The ergonomically shaped headrest ensures that patients are placed in a prone position in a safe and pressure relieving manner and that the tracheal tube remains kink-free during surgery. The headrest mirror enables the patient to be monitored by means of the high openings on the side of the headrest.

Patient safety


The individually and hygienically packed (non-sterile) headrest is only intended for single use with adults. This helps prevent crosscontamination.

Accessibility

The patient's head and face as well as the tracheal tube are all easily accessible by means of the openings on the sides.

Patient comfort

The ergonomic headrest made out of rigid PUR foam reduces the pressure on the patient's head and neck.







Single-use



Prone position headrest with mirror

For pressure-relieving positioning of patients









Features

- □ Lightweight due to the polyurethane foam (PUR), and at the same time, a high degree of shape stability
- abla The ergonomic design assists in stabilising the head without pressure building up, even during long medical surgeries
- $\ensuremath{\,\smile\,}$ The patient's head and neck are given proper support and the pressure on the patient's face is minimised
- $\ensuremath{\,\widehat{\,}}$ The openings on the sides make it possible for the endotracheal tube to remain free of
- $\ensuremath{\,\smile\,}$ Patients can be easily monitored using the headrest mirror, which comes included
- $\begin{tabular}{ll} ∇ & For use with adults \\ \end{tabular}$
- if the headrest has not been contaminated or damaged during temporary storage
- igcup Free of latex and PVC

Version	Dimension W x H x T	Colour	Pieces per box	Code	Product Code
With mirror	24 x 15 x 30 cm	Creamy-white	8	EAN 4250478428592	N14200

